

Development of Analysis Techniques for Dynamic Magnetic Resonance Imaging of the Knee Application to Tibiofemoral Kinematics



Motivation: Relevance of knee kinematics

- Knee osteoarthritis (OA) affects 10% men and 13% of women aged 60 and older (Zhang et al., 2010)
- Altered tibiofemoral kinematics in OA patients can potentially accelerate disease progression (Farrokhi et al., 2014)
- Even in ACL-deficient knees without instability symptoms, knee kinematics are altered (Yang et al., 2018)
- Understanding tibiofemoral kinematics is key to assessing knee joint function






Tibiofemoral kinematics

Six Degrees of freedom 3 Rotations

- Flexion-Extension
- Medial-Lateral
- Abduction-Adduction

3 Translations


- Anterior-Posterior
- Meidal-Lateral
- Compression-Distraction




(Credit: Gray's Anatomy 42nd ed.)

Universitätsklinikum jena

The Setup









Knee flexion-extension device

- Repeated active open-chain knee flexion-extension exercise
- 5 volunteers
- Guided by the beat of a metronome at 60 bpm
- Once without additional weight and once with 12 kgs
- Each scan lasted a duration of 160 seconds
- Range of motion achieved within the scanner ~(30-46)°



IULIUS WOLFF

BIH Center for Musculoskeletal


Biomechanics and Regeneration




Image sequence and acquisition

universitätsklinikumjena

Reconstructed frames






universitätsklinikumjena

Methods: Edge Detection

Hystereses Thresholding (Pixels and their gradient magnitude)




Low Threshold: lower values are discarded (not an edge)
High Threshold: values higher are 'Strong Edges'
In Between: 'Weak Edges' Only kept if connected to a strong edge

Strength of gaussian blur






Methods: Edge Detection





Methods: Connected Component Labelling



Treated the series
 of 2d arrays as a
 3d volume.

• 3 × 3 × 3 structuring element

• 26-way connectivity




Skeletonization

+

Remove small objects




Methods: Example sets






Methods: Issues with raw coordinates



Universitätsklinikumjena

SET R \Longrightarrow SET P




Methods: Sorting and Downsampling

- 1. Sorting:
- most distal and posterior point for the tibia, most proximal for the femur
- nearest neighbor approach to create an ordered path
- 2. Cubic Spline Interpolation:
- create a smooth, continuous curve that passes through all points in ordered set R
- parameterize the curve using arc length
- 3. Downsampling and point placement
- Equal interval = total arc length of spline / desired number of points
- place points along the spline at these equal arc length intervals




Methods: Alignment of Set P to Q





Methods: Rigid Transformation for Edge Tracking



1. A point in set P(x,y) undergoes the following transformation:

$$\begin{bmatrix} x' \\ y' \end{bmatrix} = \mathbf{R}(\phi) \begin{bmatrix} x \\ y \end{bmatrix} + \begin{bmatrix} \Delta x \\ \Delta y \end{bmatrix}$$

where, R(ϕ) is the 2D rotation matrix: $\begin{bmatrix} \cos \phi & -\sin \phi \\ \sin \phi & \cos \phi \end{bmatrix}$


And $(\Delta x, \Delta y)$ are the translations in the sagittal plane.

- Optimize for: $(\Delta x, \Delta y, \phi)$
- 2. A cost function that quantifies the misalignment between each (x', y') and its nearest neighbour in set Q.

$$C(\Delta x, \Delta y, \phi) = \sum_{p=1}^{N} \min_{q \in Q} \left(\sqrt{\left(x_q - x'_p\right)^2 + \left(y_q - y'_p\right)^2} \right)$$




Methods: Rigid Transformation for Edge Tracking




- Frame-by-frame transformation. Each transformed frame is used as reference for the next.
- Nelder-Mead simplex algorithm ("fmin search")
- Parameters:
- Initial guess: [0,0,0]
- Maximum iterations: 1000
- Function tolerance and parameter tolerance: 1e-8
- Once the set of parameters are obtained, any shape drawn on one frame on the bone will automatically be transformed to the rest of the frames




Methods: Applying transformation







Parameter extraction and tracking




Distance measurement



Angle measurement





Method: Long axis using PCA





Results









Thank you for your attention!





Supplementary materials

