

# The world of dynamic knee MRI

Different ways in which kinematic  
parameters are extracted

Jena | 00.00.2024

Aayush Nepal

## List of papers:

- Oeveleen , 2024 -> Requires one full lowerlimb scan.
  - Measured: 3D meniscal position and deformation at 0°, 30°, 60°, 90° flexion
  - Method: Static 3D MRI scans at each angle + full lowerlimb scan
  - Used statistical shape models to predict meniscal geometry
  - Calculated point-to-point distances between flexion angles
- Cance , 2024 ->
  - Measured: Direction and magnitude of discoid lateral meniscus shift
  - Method: Two static MRI scans - after 'clock' and 'pop' movements
  - Compared meniscus position between scans to classify instability type
  - Validated findings with arthroscopic examination


# The relation between meniscal dynamics and tibiofemoral kinematics

A. Van Oevelen<sup>1,2,3</sup>, M. Peiffer<sup>1,2</sup>, A. Chevalier<sup>4</sup>, J. Victor<sup>1,2</sup>, G. Steenackers<sup>3</sup>,  
E. Audenaert<sup>1,2,3,5</sup> & K. Duquesne<sup>1,2,6</sup>

- Proton density weighted 3D sequence
- (0.3571, 0.3571, 1.5 ) mm
- Additional coronal overview scans
- (0.7308, 0.7308, 3)
- This overview scan needs the full lowerlimb to create full bone models later

- Medial meniscus anterior horn, lateral meniscus anterior horn and lateral meniscus posterior horn displacement at different degrees of knee flexion (0, 30, 60, 90)

- Subjects positions in lateral decubitus (lying on their sides) on the MRI table
- Foot of the dominant leg pressed against a wooden bar
- A weight of 25 kgs attached to the wooden bar to simulate knee joint loading
- Scans were taken at four angles.
- Knee flexion angle was approximated using a goniometer\*\*
- Additional coronal scan taken at 0 degree flexion



- Use of Statistical Shape Models
- Use the overview scan and SSM .. Fit them to the MRI data to obtain full bone models.
- Use the full bone models to analyse the partial scans.

- Proton density weighted 3D sequence
- (0.3571, 0.3571, 1.5 ) mm
- Additional coronal overview scans
- (0.7308, 0.7308, 3)

- Subjects positions in lateral decubitus (lying on their sides) on the MRI table
- Foot of the dominant leg pressed against a wooden bar
- A weight of 25 kgs attached to the wooden bar to simulate knee joint loading
- Scans were taken at four angles.
- Knee flexion angle was approximated using a goniometer\*\*
- Additional coronal scan taken at 0 degree flexion







